


## Certified RNase, DNase, and Pyrogen Free Tips

For use with:	<b>Eppendorf, Pipetman, Gilson/Rainin, Finnpipette, Socorex, Titertek Digital, Rainin EDP2, Costar Onepette 8 and 12 and Volac Pipettors</b>					
<b>C110-96RNS</b>	1-200µl	10 Racks of 96, Sterile Natural	47.39	42.25	39.07	36.53




110-N

For use with:	<b>Eppendorf, Pipetman, Gilson/Rainin, Finnpipette, SMI, Socorex, Nichiryo, Costar, SMI Airpettor, and Volac Pipettors</b>					
<b>C111-R100NS</b>	101-1000µl	10 Racks of 100, Sterile Natural	51.94	45.83	43.49	41.28



111-N

For use with:	<b>Eppendorf, Nichiryo, Oxford-Benchmate, Socorex, and Titertek Plus</b>					
<b>C114-R100S</b>	0.5-10µl	10 Racks of 100, Sterile	47.39	42.25	39.07	36.53



114


# Filter Tips for PCR

## Certified RNase / DNase, Pyrogen, and Pyrogen Free

See also Rose F-Tips

Number	Description	Price per unit, in quantities of:	
		1-9	10 plus

For use with but not limited to:	<b>Oxford, Pipetman, Socorex Pipettors. Macro tip with small micro hydrophobic filter offers secure protection for your pipettor from corrosive aerosols. Macro Tip</b>			
<b>F090-100S-5000</b>	5000 $\mu$ l	Bag of 100, Sterile, Natural	22.75	18.75
<b>FA090-5000</b>	5000 $\mu$ l	Macro Tip, Bag of 100, Autoclavable, Natural	23.40	18.85



F090-5000

For use with but not limited to:	<b>Ultra Micro "Nano Tip" for use with Ultra Micro Style Pipettors including Eppendorf, Biohit Proline, Jencons Sealpette, Pipetman, Nichiryo, Oxford Benchmate</b>			
<b>F101-96RS-5</b>	5 $\mu$ l	10 Racks of 96, Sterile	86.91	70.85
<b>F101-S-5</b>	5 $\mu$ l	Single Wrapped, 10 pkg/100	120.30	96.50
<b>FA101-96R-5</b>	5 $\mu$ l	10 Racks of 96, Autoclavable	74.75	63.70
<b>FA101-5</b>	5 $\mu$ l	Bag of 1000, Autoclavable	55.58	46.31



2.5  $\mu$ l Reference mark

F101-5

For use with but not limited to:	<b>Ultra Micro "Nano Tip" for use with Ultra Micro Style Pipettors including Pipetman, Biohit Proline, Eppendorf, Jencons Sealpette, Labpette, Nichiryo, Oxford Benchmate, Rainin, Socorex</b>			
<b>FH102-96RS-10</b>	10 $\mu$ l	10 Racks of 96, Sterile, Hinged Lids	89.45	68.25
<b>F102-S-10</b>	10 $\mu$ l	Single Wrapped, 10 pkg/100	10.92	8.78
<b>FAH102-96R-10</b>	10 $\mu$ l	10 Racks of 96, Autoclavable, Hinged Lids	77.22	62.40
<b>F102-10</b>	10 $\mu$ l	Bag of 1000, Autoclavable	51.69	41.60



2.5  $\mu$ l Reference mark


F102-10

For use with but not limited to:	<b>Micro-capillary tip for gel loading, (O.D.= 0.61mm)</b>			
<b>Gel Loading Style</b>	<b>Pipetman P-2, 0.1-2 <math>\mu</math>l and P-10, 0.1-10 <math>\mu</math>l</b>			
	<b>Oxford Benchmate 0.5-10 <math>\mu</math>l, Biohit Proline Electronic Single Channel 0.2 - 10 <math>\mu</math>l</b>			
<b>F103-R204S-10</b>	10 $\mu$ l	Ultra Micro Gel Loading Tip, 5 Racks of 204, Sterile	100.75	90.35
<b>F103-S-10</b>	10 $\mu$ l	Ultra Micro Gel Loading Tip, Single Wrapped, Bag of 100	17.23	12.68



F103-10

Number	Description	Price per unit, in quantities of:	
		1-9	10 plus
For use with but not limited to:	<b>Biohit Proline: Electronic Single Channel 2-10<math>\mu</math>l, Oxford Benchmate 0.5-10<math>\mu</math>l Pipetman P-2 0.1-2<math>\mu</math>l and P-10 0.5-10<math>\mu</math>l, Jencons Sealpette</b>		
<b>F104-96RS-10</b>	10 $\mu$ l Ultra Micro Tip, 2 $\mu$ l Reference Mark, 10 Racks of 96, Sterile	86.91	70.85
<b>FA104-96R-10</b>	10 $\mu$ l Ultra Micro Tip, 10 Racks of 96, 2 $\mu$ l Reference Mark, Autoclavable	74.75	63.70
<b>FA104-10</b>	10 $\mu$ l Ultra Micro Tip, Bag of 1000, 2 $\mu$ l Reference Mark, Autoclavable, Bag of 1000	55.58	46.31




**F104-10**

For use with but not limited to:	<b>Biohit Proline 0.5-10<math>\mu</math>l, Eppendorf, Finnpiquette Digital 0.5-10<math>\mu</math>l Titertek/Labsystems Multi-channel 0.5-10<math>\mu</math>l, Jencons Sealpette</b>		
<b>F106-96RS-10</b>	10 $\mu$ l Filter Tip, 10 Racks of 96, Sterile	86.91	70.85
<b>F106-S-10</b>	10 $\mu$ l Filter Tip, Single Wrapped, Sterile, 10 x 100/bag	120.30	96.50
<b>FA106-96R-10</b>	10 $\mu$ l Filter Tip, 10 Racks of 96, Autoclavable	74.75	63.70
<b>FA106-10</b>	10 $\mu$ l Filter Tip, Bag of 1000, Autoclavable	55.58	46.31



**F106-10**

For use with but not limited to:	<b>Biohit Proline 0.5-10<math>\mu</math>l, Finnpiquette Digital 0.5-10<math>\mu</math>l, Titertek/Labsystems Multi-channel 0.5-10<math>\mu</math>l, Jencons Sealpette</b>		
<b>F108-96RS-10</b>	10 $\mu$ l Filter Tip, 10 Racks of 96, Sterile	86.91	70.85
<b>F108-S-10</b>	10 $\mu$ l Filter Tip, Single Wrapped, Sterile, 10 x 1000/bag	120.30	96.50
<b>FA108-96R-10</b>	10 $\mu$ l Filter Tip, 10 Racks of 96, Autoclavable	74.75	63.70
<b>FA108-10</b>	10 $\mu$ l Filter Tip, Bag of 1000, Autoclavable	55.58	46.31



**F108-10**

For use with but not limited to:	<b>Biohit Proline, Costar Onepette, Oxford Benchmate 10-50<math>\mu</math>l and Multi-channel 5-50<math>\mu</math>l Titertek/Labsystems Multi-channel 5-50<math>\mu</math>l, Jencons Sealpette</b>		
<b>F108-96RS-100</b>	100 $\mu$ l Filter Tip, 10 Racks of 96, Sterile	86.91	70.85
<b>F108-S-100</b>	100 $\mu$ l Filter Tip, Single Wrapped, Sterile, 10 x 100/bag	120.30	96.50
<b>FA108-96R-100</b>	100 $\mu$ l Filter Tip, 10 Racks of 96, Autoclavable	74.75	63.70
<b>FA108-100</b>	100 $\mu$ l Filter Tip, Bag of 1000, Autoclavable	55.58	46.31



**F108-100**

Number	Description	Price per unit, in quantities of:	
		1-9	10 plus
For use with but not limited to:	<b>Finnpipette 40-200<math>\mu</math>l, Biohit Proline, Oxford Benchmate 40-200<math>\mu</math>l and Multi-channel 40-200<math>\mu</math>l, Pipetman P-200 20-200<math>\mu</math>l Titertek/Labsystems Multi-channel 50-300<math>\mu</math>l</b>		
<b>F108-96RS-200</b>	200 $\mu$ l Filter Tip, 10 Racks of 96, Sterile	86.91	70.85
<b>F108-S-200</b>	200 $\mu$ l Filter Tip, Single Wrapped, Sterile, 10 x 100/bag	120.31	96.50
<b>FA108-96R-200</b>	200 $\mu$ l Filter Tip, 10 Racks of 96, Autoclavable	74.75	63.70
<b>FA108-200</b>	200 $\mu$ l Filter Tip, Bag of 1000, Autoclavable	55.58	46.31




F108-200

For use with but not limited to:	<b>Biohit Proline: Electronic Single Channel 50-200<math>\mu</math>l, Manual Single Channel 50-200<math>\mu</math>l, Finnpipette 40-200<math>\mu</math>l, 50-300<math>\mu</math>l, Oxford Benchmate 40-200<math>\mu</math>l Pipetman P-200 50-200<math>\mu</math>l Socorex (822, 852) 20-200<math>\mu</math>l, Socorex (821, 851, 853) 50-200<math>\mu</math>l Titertek 50-200<math>\mu</math>l, 50-300<math>\mu</math>l, Eppendorf Multi-channel 0 - 300<math>\mu</math>l</b>		
<b>F109-96RS-300</b>	300 $\mu$ l Filter Tip, 10 Racks of 96, Sterile	88.21	74.75
<b>F109-S-300</b>	300 $\mu$ l Filter Tip, Single Wrapped, 10 x 100/bag	120.30	96.50
<b>FA109-96R-300</b>	300 $\mu$ l Filter Tip, 10 Racks of 96, Autoclavable	74.75	63.70
<b>FA109-300</b>	300 $\mu$ l Filter Tip, Bag of 1000, Autoclavable	55.58	46.31




F109-300

For use with but not limited to: <i>Extra Long Reaching, Extended Length Tip</i>	<b>Biohit Proline, Costar, Eppendorf Digital 10-100<math>\mu</math>l Finnpipette 5-40<math>\mu</math>l, Jencons Sealpette, MLA, Nichiryo, Oxford Benchmate Pipetman P-100 10-100<math>\mu</math>l &amp; P200 50-200<math>\mu</math>l Socorex 811 &amp; 822 10-100<math>\mu</math>l, Socorex 821 &amp; 822 20-200<math>\mu</math>l</b>		
<b>F110-NXL-R204S-200</b>	200 $\mu$ l Extra Long Tip, 8 Racks of 204, Sterile, Natural	146.71	126.26
<b>F110-NXL-S-200</b>	200 $\mu$ l Extra Long Tip, Single Wrapped, Sterile, 10 x 100/Bag, Natural	129.70	108.60
<b>FA110-NXL-R204-200</b>	200 $\mu$ l Extra Long Tip, 8 Racks of 204, Autoclavable	127.21	113.10
<b>FA110-NXL-200</b>	200 $\mu$ l Extra Long Tip, Bag of 1000, Autoclavable	57.72	50.70



F110-NXL-200

For use with but not limited to:	<b>Biohit Proline: Electronic Single Channel 50-1000<math>\mu</math>l Manual Single Channel 200-1000<math>\mu</math>l, Eppendorf 100-1000<math>\mu</math>l Pipetman P-1000 100-1000<math>\mu</math>l, Socorex 811-1000, 821 &amp; 822 200-1000<math>\mu</math>l Titertek, Jencons Sealpette, Nichiryo, Oxford Benchmate</b>		
<b>F111-R100S-1000</b>	1000 $\mu$ l Extra Long Tip, 10 Racks of 100, Sterile, Natural	86.91	75.40
<b>FA111-R100-1000</b>	1000 $\mu$ l Extra Long Tip, 10 Racks of 100, Autoclavable, Natural	79.92	68.61
<b>FA111-1000</b>	1000 $\mu$ l Extra Long Tip, Bag of 1000, Autoclavable, Natural	59.15	50.70



F111-NXL-1000

Number	Description	Price per unit, in quantities of:	
		1-9	10 plus

For use with but not limited to:	<b>Eppendorf Digital 2-20<math>\mu</math>l, MLA, Costar, Pipetman P-20 2-20<math>\mu</math>l Socorex 811-10, 811-20, 822, 851, 853</b>			
<b>F113-96RNS-20</b>	20 $\mu$ l	Filter Tip, 10 Racks of 96, Sterile, Natural	86.91	70.85
<b>F113-NS-20</b>	20 $\mu$ l	Filter Tip, Single Wrapped, Sterile, 10 x 100/Bag, Natural	120.03	96.50
<b>FA113-96RN-20</b>	20 $\mu$ l	Filter Tip, 10 Racks of 96, Autoclavable, Natural	74.75	63.70
<b>FA113-N-20</b>	20 $\mu$ l	Filter Tip, Bag of 1000, Autoclavable, Natural	55.58	46.31



F113-N-20

For use with but not limited to:	<b>Biohit Proline: Electronic Single Channel 5-40<math>\mu</math>l, Finnpipette Digital 5-50<math>\mu</math>l, Jencons Sealpette, Oxford Benchmate 10-50<math>\mu</math>l, Titertek 5-40<math>\mu</math>l &amp; 5-50<math>\mu</math>l</b>			
<b>F113-96RNS-50</b>	50 $\mu$ l	Filter Tip, 10 Racks of 96, Sterile, Natural	86.91	70.85
<b>F113-NS-50</b>	50 $\mu$ l	Filter Tip, Single Wrapped, Sterile, 10 x 100/Bag, Natural	120.30	96.50
<b>FA113-96RN-50</b>	50 $\mu$ l	Filter Tip, 10 Racks of 96, Autoclavable, Natural	74.75	63.70



F113-N-50

For use with but not limited to:	<b>Eppendorf 10-100<math>\mu</math>l, Finnpipette 5-50<math>\mu</math>l, Jencons Sealpette, MLA Oxford Benchmate 40-200<math>\mu</math>l, Pipetman P-100 10-100<math>\mu</math>l</b>			
<b>F113-96RNS-100</b>	100 $\mu$ l	Filter Tip, 10 Racks of 96, Sterile, Natural	86.91	70.85
<b>F113-NS-100</b>	100 $\mu$ l	Filter Tip, Single Wrapped, Sterile, 10 x 100/Bag, Natural	120.30	96.50
<b>FA113-96RN-100</b>	100 $\mu$ l	Filter Tip, 10 Racks of 96, Autoclavable, Natural	74.75	63.70
<b>FA113-N-100</b>	100 $\mu$ l	Filter Tip, Bag of 1000, Autoclavable, Natural	55.58	46.31



F113-N-100


For use with but not limited to:	<b>Biohit Proline: Manual Single Channel 50-200<math>\mu</math>l Electronic Single Channel 10-250<math>\mu</math>l Finnpipette 40-200<math>\mu</math>l, Finnpipette 50-300<math>\mu</math>l, Nichiryo 0.5 - 10<math>\mu</math>l Pipetman P200 50-200<math>\mu</math>l, Titertek 40-200<math>\mu</math>l &amp; 50-300<math>\mu</math>l</b>			
<b>F113-96RNS-150</b>	150 $\mu$ l	Filter Tip, 10 Racks of 96, Sterile, Natural	86.91	70.85
<b>F113-NS-150</b>	150 $\mu$ l	Filter Tip, Single Wrapped, 10 x 100/Bag, Natural	120.30	96.50
<b>FA113-96RN-150</b>	150 $\mu$ l	Filter Tip, 10 Racks of 96, Autoclavable, Natural	74.75	63.70



F113-N-150

Number	Description	Price per unit, in quantities of:	
		1-9	10 plus

For use with but not limited to:	<b>Eppendorf 0.5-10<math>\mu</math>l, Oxford Benchmate 0.5-10<math>\mu</math>l, Socorex 822-0002, 0.4-2<math>\mu</math>l</b>			
<b>F114-R100S-10</b>	10 $\mu$ l	Ultra Micro Tip, 10 Racks of 100, Sterile, Natural	86.91	70.85
<b>F114-S-10</b>	10 $\mu$ l	Ultra Micro Tip, Single Wrapped, Sterile, 10 x 100/Bag	120.03	9.65
<b>FA114-R100-10</b>	10 $\mu$ l	Ultra Micro Tip, 10 Racks of 100, Autoclavable	74.75	63.70
<b>FA114-10</b>	10 $\mu$ l	Ultra Micro Tip, Bag of 1000, Autoclavable	55.58	46.31



**F114-10**

For use with but not limited to:	<b>Oxford Benchmate 0.5-10<math>\mu</math>l Pipetman P20 2-20<math>\mu</math>l (P20 Requires adapter, included with tips).</b>			
<b>F114-R100S-20</b>	20 $\mu$ l	Ultra Micro Tip, 10 Racks of 100, Sterile	86.91	70.85
<b>F114-S-20</b>	20 $\mu$ l	Ultra Micro Tip, Single Wrapped, Sterile, 10 x 100/bag	120.30	96.50
<b>FA114-R100-20</b>	20 $\mu$ l	Ultra Micro Tip, 10 Racks of 100, Autoclavable	74.75	63.70
<b>FA114-20</b>	20 $\mu$ l	Ultra Micro Tip, Bag of 1000, Autoclavable	55.58	46.31




**F114-20**

For use with but not limited to:	<b>Biohit Proline, Eppendorf 10-100<math>\mu</math>l Model 4850 Electronic &amp; Model 5221 EDOS. Model 4810 adjustable volume, Model 4800 fixed volume. Model 3190 standard fixed volume, Model 4703 tri-volume. Jencons Sealpette, MLA, Model 2000, Model 2100</b>			
<b>F115-96RNS-100</b>	100 $\mu$ l	Ultra Micro Tip, 10 Racks of 96, Sterilized, Natural	86.91	70.85
<b>F115-NS-100</b>	100 $\mu$ l	Ultra Micro Tip, Single Wrapped, Sterile, 10 x 100/Bag, Natural	120.30	96.50
<b>FA115-96RN-100</b>	100 $\mu$ l	Ultra Micro Tip, 10 Racks of 96, Autoclavable, Natural	74.75	63.70
<b>FA115-N-100</b>	100 $\mu$ l	Ultra Micro Tip, Bag of 1000, Autoclavable, Natural	55.58	46.31



**F115-N-100**

For use with but not limited to:	<b>Eppendorf 10-100<math>\mu</math>l, Oxford Benchmate 40-200<math>\mu</math>l, Finnpiquette, MLA Extra-long reaching tips. Ideal for pipetting from test tubes.</b>			
<b>F115-NXL-R100S-200</b>	200 $\mu$ l	Extra Long Tip, 8 Racks of 100, Sterile, Natural	80.41	68.25
<b>F115-NXL-S-200</b>	200 $\mu$ l	Extra Long Tip, Single Wrapped, Sterile, 10 x 100/Bag, Natural	120.30	96.50
<b>FA115-NXL-R100-200</b>	200 $\mu$ l	Extra Long Tip, 8 Racks of 100, Autoclavable, Natural	66.95	57.20
<b>FA115-NXL-200</b>	200 $\mu$ l	Extra Long Tip, Bag of 1000, Autoclavable, Natural	61.91	52.00



**F115-NXL-200**


Number	Description	Price per unit, in quantities of:	
		1-9	10 plus

For use with but not limited to:	<b>Eppendorf 350<math>\mu</math>l, Model 4850 – Electronic, Model 5221 - EDOS</b>			
<b>F116-R100NS-350</b>	350 $\mu$ l	Filter Tip, 10 Racks of 100, Sterile, Natural	86.91	75.40
<b>F116-NS-350</b>	350 $\mu$ l	Filter Tip, Single Wrapped, Sterile, 10 x 100/Bag, Natural	124.60	106.00
<b>FA116-R100N-350</b>	350 $\mu$ l	Filter Tip, 10 Racks of 100, Autoclavable, Natural	79.92	68.61
<b>FA116-N-350</b>	350 $\mu$ l	Filter Tip, Bag of 1000, Autoclavable, Natural	59.15	50.70




**F116-N-350**

For use with but not limited to:	<b>Eppendorf 500<math>\mu</math>l, Biohit Proline, Jencons Sealpette</b>			
<b>F116-NXL-R100S-500</b>	500 $\mu$ l	Filter Tip, 10 Racks of 100, Sterile, Natural	86.91	75.40
<b>F116-NXL-S-500</b>	500 $\mu$ l	Filter Tip, Single Wrapped, Sterile, 10 x 100/bag, Natural	124.60	106.00



**F116-NXL-500**

For use with but not limited to:	<b>Eppendorf, Biohit Proline, Finnpiquette, Jencons Sealpette, Nichiryo, Oxford Benchmate, Socorex, Titertek</b>			
<b>F116-NXL-R100S-1000</b>	1000 $\mu$ l	Filter Tip, 10 Racks of 100, Sterile, Natural	86.91	75.40
<b>F116-NXL-S-1000</b>	1000 $\mu$ l	Filter Tip, Single Wrapped, 10 x 100/Bag, Natural	124.60	106.00




**F116-NXL-1000**

For use with but not limited to:	<b>Biohit Proline - Manual Single Channel 5-50<math>\mu</math>l Biohit Proline - Electronic Single Channel 5-100<math>\mu</math>l Biohit Proline - Manual Multi-channel 5-100<math>\mu</math>l</b>			
<b>F117-96RNS-100</b>	100 $\mu$ l	Filter Tip, 10 Racks of 100, Sterile, Natural	86.91	70.85
<b>F117-NS-100</b>	100 $\mu$ l	Filter Tip, Single Wrapped, 10 x 100/Bag, Natural	120.30	96.50
<b>FA117-96RN-100</b>	100 $\mu$ l	Filter Tip, 10 Racks of 100, Autoclavable, Natural	74.75	63.70
<b>FA117-N-100</b>	100 $\mu$ l	Filter Tip, Bag of 1000, Autoclavable, Natural	55.58	46.31



**F117-N-100**


For use with but not limited to:	<b>Biohit Proline - Manual Single Channel 50-200<math>\mu</math>l Biohit Proline - Electronic Single Channel 10-250<math>\mu</math>l Biohit Proline - Electronic Multi-channel 25-250<math>\mu</math>l</b>			
<b>F117-96RNS-200</b>	200 $\mu$ l	Filter Tip, 10 Racks of 96, Sterile, Natural	86.91	70.85
<b>F117-NS-200</b>	200 $\mu$ l	Filter Tip, Single Wrapped, 10 x 100/Bag, Natural	120.30	96.50
<b>FA117-96RN-200</b>	200 $\mu$ l	Filter Tip, 10 Racks of 96, Autoclavable, Natural	74.75	63.70
<b>FA117-N-200</b>	200 $\mu$ l	Filter Tip, Bag of 1000, Autoclavable, Natural	55.58	46.31



**F117-N-200**


Number	Description	Price per unit, in quantities of:	
		1-9	10 plus

For use with but not limited to:	<b>Biohit Proline - Manual Multi Channel 50-250<math>\mu</math>l, &amp; 50-300<math>\mu</math>l</b> <b>Biohit Proline - Electronic Single Channel 10-250<math>\mu</math>l</b> <b>Biohit Proline - Electronic Multi-channel 25-250<math>\mu</math>l</b> <b>Jencons Sealpette, Manual &amp; Electronic, Single &amp; Multi-channels</b>			
<b>F117-NXL-96RNS-300</b>	300 $\mu$ l	Filter Tip, 10 Racks of 96, Sterile, Natural	86.91	70.85
<b>F117-NXL-S-300</b>	300 $\mu$ l	Filter Tip, Single Wrapped, 10 x 100/Bag, Natural	120.30	96.50
<b>FA117-NXL-96R-300</b>	300 $\mu$ l	Filter Tip, 10 Racks of 96, Autoclavable, Natural	74.75	63.70
<b>FA117-NXL-300</b>	300 $\mu$ l	Filter Tip, Bag of 1000, Autoclavable, Natural	55.58	46.31




F117-NXL-300

For use with but not limited to: <b>Large Orifice Tip!</b>	<b>Biohit Proline 50-200<math>\mu</math>l</b> <b>Costar Onepette, Eppendorf, Finnpiquette 40-200<math>\mu</math>l</b> <b>Nichiryo, Oxford Benchmate, Titertek</b> <b>Jencons Sealpette, Pipetman P200 50-200<math>\mu</math>l</b>			
<b>F118-96RNS-200</b>	200 $\mu$ l	Large Orifice Tip, 10 Racks of 96, Sterile, Natural	86.91	70.85
<b>F118-NS-200</b>	200 $\mu$ l	Large Orifice Tip, Single Wrapped, 10 x 100/Bag, Natural	120.30	96.50
<b>FA118-96RN-200</b>	200 $\mu$ l	Large Orifice Tip, 10 Racks of 96, Autoclavable, Natural	74.75	63.70
<b>FA118-N-200</b>	200 $\mu$ l	Large Orifice Tip, Bag of 1000, Autoclavable, Natural	55.58	46.31




F118-N-200

For use with but not limited to:	<b>Biohit Proline: Manual Single Channel 200-1000<math>\mu</math>l</b> <b>Electronic Single Channel 50-1000<math>\mu</math>l</b> <b>Costar, Eppendorf 100-1000<math>\mu</math>l, Finnpiquette 200-1000<math>\mu</math>l, Oxford Benchmate 200-1000<math>\mu</math>l</b> <b>Nichiryo, Titertek, Pipetman P1000 100-1000<math>\mu</math>l, Socorex 200-1000<math>\mu</math>l</b>			
<b>F119-NXL-R100S-1000</b>	1000 $\mu$ l	Large Orifice Extra Long Tip, 10 Racks of 100, Sterile, Natural	86.91	75.40
<b>FA119-NXL-R100-1000</b>	1000 $\mu$ l	Large Orifice Tip, 10 Racks of 100, Autoclavable, Natural	79.92	68.61
<b>FA119-NXL-1000</b>	1000 $\mu$ l	Large Orifice Tip, Bag of 1000, Autoclavable, Natural	59.15	50.70



F119-NXL-1000

For use with but not limited to:	<b>Oxford Benchmate 200-1000<math>\mu</math>l</b> <b>Pipetman P1000 100-1000<math>\mu</math>l, Titertek</b>			
<b>F122-96RS-1250</b>	1250 $\mu$ l	Filter Tip, 10 Racks of 96, Sterile	86.91	75.40
<b>F122-S-1250</b>	1250 $\mu$ l	Filter Tip, Single Wrapped, 10 x 100/Bag	124.61	105.95
<b>FA122-96R-1250</b>	1250 $\mu$ l	Filter Tip, 10 Racks of 96, Autoclavable	79.92	68.61
<b>FA122-1250</b>	1250 $\mu$ l	Filter Tip, Bag of 1000, Autoclavable	59.15	50.70




F122-1250

Number	Description	Price per unit, in quantities of:	
		1-9	10 plus
For use with but not limited to: <b>Gel Loading Style</b>	<b>Micro-capillary tip for gel loading (O.D. = 0.61mm)</b> Eppendorf, 0.5-10 $\mu$ l, Pipettors, Nichiryo, 0.5-10 $\mu$ l, Oxford Benchmate 0.5-10 $\mu$ l Socorex 811 5 $\mu$ l, Socorex 822 0.4 $\mu$ l and 1-10 $\mu$ l		
<b>F123-S-10</b>	0.1-10 $\mu$ l	Ultra Gel Loading Tip, Single Wrapped, Sterile, 100/Bag	17.23    15.60




F123-10

For use with but not limited to:	<b>Biohit Proline Electronic Multi-channel 50-1200<math>\mu</math>l</b> <b>Jencons Sealpette Electronic Multi-channel 50-1200<math>\mu</math>l</b>		
<b>F126-96RS-1200</b>	1200 $\mu$ l	Filter Tip, 10 Racks of 96, Sterile	86.91    70.85
<b>F126-S-1200</b>	1200 $\mu$ l	Filter Tip, Single Wrapped, 10 x 100/Bag	120.25    96.53
<b>FA126-96R-1200</b>	1200 $\mu$ l	Filter Tip, 10 Racks of 96, Autoclavable	74.75    63.70
<b>FA126-1200</b>	1200 $\mu$ l	Filter Tip, Bag of 1000, Autoclavable	55.58    46.31




F126-1200

For use with but not limited to:	<b>Biohit Proline Manual Single-channel 50-200<math>\mu</math>l</b> <b>Biohit Proline Manual Multi-channel 50-250<math>\mu</math>l</b> <b>Biohit Proline Manual Multi-channel 50-300<math>\mu</math>l</b> <b>Biohit Proline Electronic Single-channel 10-250<math>\mu</math>l</b> <b>Biohit Proline Electronic Multi-channel 25-250<math>\mu</math>l</b> <b>Costar Onepette, Eppendorf, Finnpiette, Jencons Sealpette Manual &amp; Electronic Single and Multi-channel, Pipetman, Rainin, Titertek</b>		
<b>F169-96RS-300</b>	300 $\mu$ l	Filter Tip, 10 Racks of 96, Sterile, Natural	86.91    70.85
<b>F169-S-300</b>	300 $\mu$ l	Filter Tip, Single Wrapped, 10 x 100/Bag, Natural	120.25    96.20



F169-300

For use with but not limited to:	<b>Socorex 811 - 100<math>\mu</math>l, 821-200, 822-0100, 822-0200</b> <b>Socorex 852-100, 822-200 (Multi-channel)</b>		
<b>F177-96RS-100</b>	100 $\mu$ l	Filter Tip, 10 Racks of 96, Sterile, Natural	86.91    75.40
<b>FA177-96R-100</b>	100 $\mu$ l	Filter Tip, 10 Racks of 96, Autoclavable, Natural	79.92    68.61
<b>FA177-100</b>	100 $\mu$ l	Filter Tip, Bag of 1000, Autoclavable, Natural	59.15    50.70



F177-100

## High Quality Racked Filter Tips: RoseF-Tips™

These low adhesion filter tips adopt the very latest micropore technology to ensure that all aerosols and liquids are prevented from entering the pipette. This technology eliminates the need for additives such as sodium, found in self-sealing type filters, which has the disadvantage of trapping samples in the tip once the filter is sealed. **RoseF-Tip™** Filter Tips use virgin polyethylene filters without additives (which is a potential source of contamination), and the micropore technology ensures that the filters will never trap samples in the tip. Should the filter become wet due to inadvertent over-pipetting, it will act as a barrier to the passage of liquid, yet the sample in the tip may still be dispensed in the normal way, and can be re-used without fear of contamination. All **RoseF-Tip™** Filter Tips are metal free and are ideal for DNA amplification (PCR) and microbiological applications, including the handling of radioactive, bio-hazardous, and corrosive samples. **Simply the best filter tip available!**

**Certified Sterile, RNase, DNase & Pyrogen Free**



**Hydrophobic sterile aerosol barrier filter tips**

**RoseF-Tip™ Aerosol Barrier FilterTips for research, diagnostic, and routine applications.**

**Certified Sterile, RNase, DNase & Pyrogen Free.**

Catalog No.	Working Volume	Color	Package Size	Price
2-001-96-0	0.1-10µl	Natural	10 Racks (96 Tips) / 960 Tips	55.25
2-051-96-0	0.5-10µl	Natural	10 Racks (96 Tips) / 960 Tips	55.25
2-114-96-0	1-20µl	Natural	10 Racks (96 Tips) / 960 Tips	55.25
2-116-96-0	1-30µl	Natural	10 Racks (96 Tips) / 960 Tips	55.25
2-111-96-0	1-50µl	Natural	10 Racks (96 Tips) / 960 Tips	55.25
2-112-96-0	1-100µl	Natural	10 Racks (96 Tips) / 960 Tips	55.25
2-113-96-0	1-180µl	Natural	10 Racks (96 Tips) / 960 Tips	55.25
2-110-96-0	1-200µl	Natural	10 Racks (96 Tips) / 960 Tips	55.25
2-115-96-0	10-230µl	Natural	10 Racks (96 Tips) / 960 Tips	55.25
2-201-96-0	100-1000µl	Natural	10 Racks (96 Tips) / 960 Tips	55.25

**Compatibility Chart:** For use with, but *not* restricted to the following brand of pipettes.

0.1-10µl <b>RoseF-Tip™</b>	Eppendorf Reference & Research, Rose Witopet, Socorex, Nichiryo, Oxford Benchmate, Titertek
0.5-10µl <b>RoseF-Tip™</b>	Pipetman P2 & P10, Rose Witopet, Eppendorf Reference & Research, Finpipet, Biohit Proline, Nichiryo, Oxford Benchmate, Rainin EDP-2
20-230µl <b>RoseF-Tip™</b> (1-20µl, 1-30µl, 1-50µl, 1-100µl, 1-180µl, 1-200µl, 1-230µl)	Gilson, Eppendorf Reference & Research, Brand, Labsystems Finnpipet, Biohit Proline, Socorex, Nichiryo, MLA, Costar, Oxford Benchmate, Rainin EDP-2, Volac, Jencons, Vanguard, Rose Witopet
100-1000µl <b>RoseF-Tip™</b>	Gilson, Eppendorf Reference & Research, Brand, Labsystems Finnpipet, Biohit Proline, Socorex, Nichiryo, MLA, Costar, Oxford Benchmate, Rainin EDP-2, Volac, Vanguard, Rose Witopet



All filter tips are certified DNase/RNase and Pyrogen free. We boast the widest selection of filter tips. If you have a pipettor that you are unsure about, please call us. Our long list of users includes most universities, research hospitals, government agencies, and many other professional organizations.

### Other PCR Related Equipment and Supplies:

**0.2ml thin wall PCR tubes in a 96, 48, 32 & 24 well, microtiter format.** Certified RNase and DNase free, these tubes can be used with the dome style caps eliminating the need to use a mineral overlay when used in thermal cyclers with heated lids.

**PCR tubes in strips.** Designed for optimal fit in all thermal cyclers. Certified RNase and DNase free.