


## Transax Tips



The Transax Transfer system was designed with the environment in mind using 90% less plastic than conventional racks and packaging made from recycled paper. Also, by allowing you to reuse existing pipette tip racks, storage space is optimized and waste is minimized.

Transferring tips is easy with the Transax™ transferring device.

- Easy manual transfer
- Universal card fits all manufacturers racks
- Snaps securely onto custom designed Medax rack




### Transax™ System Starter Kits For Ultra-Micro Tips, Standard 0-200µl Tips and 100-1000µl Tips

Number	Description	Price
<b>RMTSK-001</b>	<b>Ultra-Micro Tips in Hinged Racks</b> 10 racks of 96/rack tips plus 10 refills for a total of 1920 tips	54.51
<b>RM-TR-001</b>	Refill for above, 10µl, 10 x 96 tips/tray	25.12
<b>RMTSK-002</b>	<b>0 – 200µl Tips in Hinged Racks</b> 10 racks of 96 tips/rack plus 10 refills for a total of 1920 tips	49.77
<b>RM-TR-002</b>	Refills for above, 200µl. 10 x 96 tips/tray	23.22
<b>RMTSK-003</b>	<b>100 – 1000µl Tips in Hinged Racks</b> 12 racks of 70 tips/rack plus 12 refills for a total of 1680 tips	62.56
<b>RM-TR-003</b>	Refills for above, 1000 µl, 10 x 70 tips/tray	29.39

## Rose Simplest™ Sliding Reload Tips System

Now you can buy your popular Rose Brand Pipette Tips to use with the simplest Tip Reloading System ever made. It is so simple that it does not require instructions. Simply take off a tray of 96 tips from your pack of 10 racks and slide it into the groove of the Simplest Holder™ and insert into the autoclavable, hinged box supplied. All this in less than 10 seconds with absolutely no fuss! It's inexpensive, re-useable, and does not require expensive accessory or hardware.



Along with the simplicity of Simplest™ comes economical benefits; you no longer have to sacrifice your budget for convenience and environmental consciousness. A case of Simplest™ Tips would typically cost you less than a case of Racked Tips from competitors.

Number	Volume	Package Type	Price per unit, in quantities of:			
			1-10	11-30	31-50	51 plus
For use with:	<b>Eppendorf, Pipetman, Gilson/Rainin, Finnpiette, Socorex, Wiptopet, etc.</b>					
<b>Q-104-R-960</b>	0.1-10 $\mu$ l	10 Racks x 96, Natural	27.98	26.59	25.25	23.99
<b>Q-110-B-R-960</b>	1-200 $\mu$ l	10 Racks x 96, Yellow	26.81	22.88	21.55	20.74
<b>Q-110-BN-R-960</b>	1-200 $\mu$ l	10 Racks x 96, Natural	26.81	22.88	21.55	20.74
<b>Q-113-G-R-960</b>	1-200 $\mu$ l	10 Racks x 96, Yellow	26.81	22.88	21.55	20.74
<b>Q-113-GN-R-960</b>	1-200 $\mu$ l	10 Racks x 96, Natural	26.81	22.88	21.55	20.74
<b>Q-9434-SET</b>	Simplest™ Sliding Reloader in Hinged Box. Pack of 10		48.75			

## High Quality Racked Filter Tips: RoseF-Tips™

RoseF-Tip™ Aerosol Barrier FilterTips for research, diagnostic, and routine applications.

Certified Sterile, RNase, DNase & Pyrogen Free.

Catalog No.	Working Volume	Color	Package Size	Price
2-001-96-0	0.1-10µl	Natural	10 Racks (96 Tips) / 960 Tips	55.25
2-051-96-0	0.5-10µl	Natural	10 Racks (96 Tips) / 960 Tips	55.25
2-114-96-0	1-20µl	Natural	10 Racks (96 Tips) / 960 Tips	55.25
2-116-96-0	1-30µl	Natural	10 Racks (96 Tips) / 960 Tips	55.25
2-111-96-0	1-50µl	Natural	10 Racks (96 Tips) / 960 Tips	55.25
2-112-96-0	1-100µl	Natural	10 Racks (96 Tips) / 960 Tips	55.25
2-113-96-0	1-180µl	Natural	10 Racks (96 Tips) / 960 Tips	55.25
2-110-96-0	1-200µl	Natural	10 Racks (96 Tips) / 960 Tips	55.25
2-115-96-0	10-230µl	Natural	10 Racks (96 Tips) / 960 Tips	55.25
2-201-96-0	100-1000µl	Natural	10 Racks (96 Tips) / 960 Tips	55.25

**Compatibility Chart:** For use with, but *not* restricted to the following brand of pipettes.

0.1-10µl <b>RoseF-Tip™</b>	Eppendorf Reference & Research, Rose Witopet, Socorex, Nichiryo, Oxford Benchmate, Titertek
0.5-10ul <b>RoseF-Tip™</b>	Pipetman P2 & P10, Rose Witopet, Eppendorf Reference & Research, Finpipet, Biohit Proline, Nichiryo, Oxford Benchmate, Rainin EDP-2
20-230ul <b>RoseF-Tip™</b> (1-20µl, 1-30µl, 1-50µl, 1-100µl, 1-180µl, 1-200µl, 1-230µl)	Gilson, Eppendorf Reference & Research, Brand, Labsystems Finnpipet, Biohit Proline, Socorex, Nichiryo, MLA, Costar, Oxford Benchmate, Rainin EDP-2, Volac, Jencons, Vanguard, Rose Witopet
100-1000µl <b>RoseF-Tip™</b>	Gilson, Eppendorf Reference & Research, Brand, Labsystems Finnpipet, Biohit Proline, Socorex, Nichiryo, MLA, Costar, Oxford Benchmate, Rainin EDP-2, Volac, Vanguard, Rose Witopet