


Rotary Evaporators – Buchi R210 / 215

Advanced and Professional Systems:

R210 and R215 with a very wide choice of Glassware, Baths, Motorized Lift System, Accessories and Automation with a choice of integrated Vacuum Pumps, Vacuum Controllers and Chillers, etc.

R215 features digital display of speed and vapor pressure/temperature. Both models are supplied with the Model B-491 Bath with digital display of temperature and having a 4 liter capacity. An optional larger 5 liter capacity bath is also available allowing for evaporation with 5 liter evaporating flask – making these systems very useful for high throughput synthesis laboratories.

Dimensions (L x H x D)	550 x 575 x 415 mm
Weight (depending on Glass configuration)	16 - 18 kg without bath
Operating voltage	100 - 240 V
Frequency	50/60 Hz
IP class	21
Power consumption	60 W (1360W with bath)
Rotation speed	20 - 280 rpm
Flask size	50 - 4000 ml
Max. flask weight	3 kg
Display: R210/R215	No/Yes
Seal of approval	UL, CSA



Standard supply includes Motor, Glassware, B-491 Heating Bath, Electronic Lift System, 2 liter Evaporation and Receiving Flask, Vapor Duct, Seal and Manual.

Ordering Information: Prices are listed in US Dollars

		Prices/USD
23022A000	R210 Advanced Rotary Evaporator with Diagonal Condenser, 4 Liter Heating Bath with Digital Display(180°C), Motorized Lift System, and all necessary glassware in standard North American joints.	5,310.75
23022A100	Same as above, with Plastic Coated Glassware	5,483.85
23022V000	R210 Advanced Rotary Evaporator with Vertical Condenser, 4 Liter Heating Bath with Digital (180°C), Motorized Lift System, and all necessary glassware in standard North American joints.	5,469.95
23022V100	Same as above, with Plastic Coated Glassware	5,643.10
23022C000	R210 Advanced Rotary Evaporator with Dry-Ice Condenser, 4 Liter Heating Bath with Digital Display (180°C), Motorized Lift System, and all necessary glassware in standard North American joints.	5,802.30
23022C100	Same as above, with Plastic Coated Glassware	5,975.45
23122V000	R215 Professional Rotary Evaporator with Diagonal Condenser, 4 Liter Heating Bath with Digital Display(180°C), Motorized Lift System, and all necessary glassware in standard North American joints.	6,251.10
23122V100	Same as above, with Plastic Coated Glassware	6,422.05
23122A000	R215 Professional Rotary Evaporator with Vertical Condenser, 4 Liter Heating Bath with Digital (180°C), Motorized Lift System, and all necessary glassware in standard North American joints.	5,959.90
23122A100	Same as above, with Plastic Coated Glassware	6,129.80
23122C000	R215 Professional Rotary Evaporator with Dry-Ice Condenser, 4 Liter Heating Bath with Digital Display (180°C), Motorized Lift System, and all necessary glassware in standard North American joints.	6,442.40
23122C100	Same as above, with Plastic Coated Glassware	6,612.30
Important Notes		
Other types of condensers are also available including double coiled and with vapor temperature options.		
Also, the 4 Liter Baths could be optionally substituted with a Larger 5 Liter Bath for use with larger evaporation flasks commonly used in custom synthesis laboratories.		

Rose Scientific Inc.
Telephone: 1-800-661-9288
sales@rosesci.com, usales@rosesci.com