

Filter Tips

Terrific introductory prices from Rose Scientific Inc.

Rose Brand Plasticware. Millions sold in North America

High Quality Racked Filter Tips:

These low adhesion filter tips adopt the very latest micropore technology to ensure that all aerosols and liquids are prevented from entering the pipette. This technology eliminates the need for additives such as sodium, found in self-sealing type filters, which has the disadvantage of trapping samples in the tip once the filter is sealed. **RoseF-Tipä** Filter Tips use virgin polyethylene filters without additives (which is a potential source of contamination), and the micropore technology ensures that the filters will never trap samples in the tip. Should the filter become wet due to inadvertent over-pipetting, it will act as a barrier to the passage of liquid, yet the sample in the tip may still be dispensed in the normal way, and can be re-used without fear of contamination. All **RoseF-Tipä** Filter Tips are metal free and are ideal for DNA amplification (PCR) and microbiological applications including the handling of radioactive, bio-hazardous and corrosive samples. **Simply the best filter tip available!**

Certified Sterile, RNase, DNase & Pyrogen Free



Hydrophobic sterile aerosol barrier filter tips

RoseF-Tipä Aerosol Barrier FilterTips for research, diagnostic and routine applications. Certified Sterile, RNase, DNase & Pyrogen Free.

Catalog No.	Working Volume	Color	Package Size	Price
2-001-96-0	0.1-10 µl	Natural	10 Racks (96 tips) / 960 Tips	55.25
2-051-96-0	0.5-10 µl (crystal)	Natural	10 Racks (96 tips) / 960 Tips	55.25
2-114-96-0	1-20 µl	Natural	10 Racks (96 tips) / 960 Tips	55.25
2-116-96-0	1-30 µl	Natural	10 Racks (96 tips) / 960 Tips	55.25
2-111-96-0	1-50 µl	Natural	10 Racks (96 tips) / 960 Tips	55.25
2-112-96-0	1-100 µl	Natural	10 Racks (96 tips) / 960 Tips	55.25
2-113-96-0	1-180 µl	Natural	10 Racks (96 tips) / 960 Tips	55.25
2-110-96-0	1-200 µl	Natural	10 Racks (96 tips) / 960 Tips	55.25
2-115-96-0	10-230 µl	Natural	10 Racks (96 tips) / 960 Tips	55.25
2-201-96-0	100-1000 µl	Natural	10 Racks (96 tips) / 960 Tips	55.25

Purchase any combination of 10 packs and receive a further 10% discount off the above special prices.

Compatibility Chart: For use with, but **not** restricted to the following brand of pipettes.

0.1-10µl RoseF-Tipä	Eppendorf Reference & Research, Rose Witopet, Socorex, Nichiryo, Oxford Benchmate, Titertek
0.5-10ul RoseF-Tipä	Pipetman P2 & P10, Rose Witopet, Eppendorf Reference & Research, Finpipet, Biohit Proline, Nichiryo, Oxford Benchmate, Rainin EDP -2
20-230ul RoseF-Tipä (1-20µl, 1-30µl, 1-50µl, 1-100µl, 1-180µl, 1-200µl, 1-230µl)	Gilson, Eppendorf Reference & Research, Brand, Labsystems Finnpipet, Biohit Proline, Socorex, Nichiryo, MLA, Costar, Oxford Benchmate, Rainin EDP -2, Volac, Jencons, Vanguard, Rose Witopet
100-1000µl RoseF-Tipä	Gilson, Eppendorf Reference & Research, Brand, Labsystems Finnpipet, Biohit Proline, Socorex, Nichiryo, MLA, Costar, Oxford Benchmate, Rainin EDP -2, Volac, Vanguard, Rose Witopet

Rose Scientific Inc.

Phone: 1-800-661-9288

Email: USAsales@rosesci.com

Fax: 1-513-942-4077

www.rosescientific.com