

# MAGNETIC STIRRER MS-01



**MAGNETIC STIRRER MS01** is a powerful and simple to use four position stainless still magnetic stirrer.

#### PRODUCT FEATURES

- Four positions for effective mixing in different size flat bottom laboratory flasks.
- Compatibility with ELMI thermostat TW-2.03
- Stainless still highly resistible body.
- Can be easily integrated in different lab processes.

#### PRODUCT BENEFITS

It is specially designed to operate with ELMI water thermostat TW-2 series. Together with TW-2.03 heating functions it becomes possible to maintain highly accurate temperature maintenance as well as mixing efficiency control. Reliability, simple operation and noiseless operation make MS01 perfect assistant in any laboratory.

## MAGNETIC STIRRER MS01 SPECIFICATIONS

Speed range	RPM	80 - 1 500
Liquid volume range for one mixing place	L	from 0.2 up to 2L
Power adapter output	$\ominus \text{---} \oplus$	DC in 12V, 1 000 mA
Dimensions (length x width x height)	mm	210 x 270 x 35
Weight	kg	1.8

### Rose Scientific Ltd.

enhancing the sciences...

4027 - 97 Street  
Edmonton, Alberta, Canada. T6E 5Y5  
Tel: 780-438-5110, 800-661-9289  
Fax: 780-462-5776, 800-570-6067

Mississauga Office:  
6581 Kitimat Road, Unit #7  
Mississauga, Ontario, Canada. L5N 3T5  
Tel: 877-826-6706  
Fax: 905-814-1978  
sales@rosesci.com