


ELISA MICROPLATE THERMOSTAT MT-2




x2



MT-2 is great accuracy thermostat for 2x different micro-plates.

Advanced double heating system which insures heating from the lead and from the platform together with the hermetic aluminum container and digital platinum sensors allow to achieve incredible heating accuracy.



SPECIFICATIONS

Time controller:	min	1-999/hold
Environmental temperature:	°C	from +4 to + 55
Relative moisture of atmosphere (at the temperature of 20 °C):	%	up to 80
Range of the temperature:		ambient- 60
Temperature of stabilization at room temperature:	°C	5
Accuracy of the temperature on the platform:	°C	0,1
Power supply	⊖ ⊕ ⊖	220V, 50Hz
Size of the device (length x width x height)	mm	310 x 270 x 115
Weight	kg	6,4

Rose Scientific Ltd.
enhancing the sciences...

4027 - 97 Street
Edmonton, Alberta, Canada. T6E 5Y5
Tel: 780-438-5110, 800-661-9289
Fax: 780-462-5776, 800-570-6067
Email: sales@rosesci.com

Mississauga Office:
6581 Kitimat Road, Unit #7
Mississauga, Ontario, Canada. L5N 3T5
Tel: 877-826-6706
Fax: 905-814-1978