


VORTEX MIXER V-3

2 TYPES OF EXCHANGEABLE VORTEX ADAPTERS

CONICAL-TYPE ADAPTER

GRID-TYPE ADAPTER



ANY KIND OF TEST TUBES
from 0.5 ml up to 50 ml



VORTEX MIXER V-3

VORTEX V-3 is a pretty looking and very easy handled vortex mixer for any type of test tubes of up to 50 ml volume.

PRODUCT FEATURES

- ▶ Touch or hands free mode.
- ▶ Easy vortex intensity regulation.
- ▶ 2 exchangeable vortexing adapters.
- ▶ Quiet and maintenance free mechanism.
- ▶ Half-sphere shape for maximal stability on the desktop.

PRODUCT BENEFITS

V-3 is simple and handy solution for current laboratory tasks.

VORTEX MIXER V-3 SPECIFICATIONS

Socket oscillation speed	RPM	50 to 3400
Socket oscillation amplitude	mm	4
Power adapter output	⊖ ⊕	DC in 12V, 500 mA
Power consumption	W	5
Dimensions (diameter x height)	mm	130 x 60
Weight	kg	0.6