

# Dry Block Heater



This block heater is a convenient answer to the laboratory's need for a small constant temperature heater that can accommodate a variety of different sizes. Heater opening is a standard SBS format, which allows heating of standard well plates (low or high profile). Our block heater features a digital control with easy set-up and simple input commands; makes this control very "user friendly". This control is single-loop and auto-tuning that automatically set PID control parameters for optimum system performance. Maintains temperature from slightly above ambient to 125°C +/- 1°C.

### Specifications

Electrical: On/Off power switch, fuse protection, 3-wire line cord with molded plug.

Accuracy: 0.1% of span.

Ambient Temperature Range: 32-149°F.

Display: Dual, 4-digit LED.

Set Point: Up/down keys. °F or °C display, user selectable.

Finish: Black

**Note: Aluminum Blocks sold separately**



Catalog Number	Watts	VAC	Dimensions (w x h x d)	Shpg wt lb (kg)	Price
099A HBC12	150	120	5-3/4 x 4-1/2 x 6-3/4	5 (2.3)	\$670.00

Catalog Number	Block Description	Shpg wt lb (kg)	Price
099A R34 562	Holds 24 (12mm x 32mm) vials	2 (0.9)	\$95.00
099A R34 561	Holds 24 (15mm x 45mm) vials	2 (0.9)	\$95.00
099A R34 563	Holds 20 (15mm or 16mm) tubes	2 (0.9)	\$95.00
099A R34 516	Holds 8 EPA vials	2 (0.9)	\$95.00

**REPEATABLE.**

**RELIABLE.**

**RESULTS.**

**Glas-Col, LLC**

711 Hulman Street., P.O. Box 2128, Terre Haute, IN 47802-0128 USA

Phone. (812) 235-6167, Fax. (812) 234-6975

E-mail: [pinnacle@glascol.com](mailto:pinnacle@glascol.com)

[www.glascol.com](http://www.glascol.com)



**Rose Scientific Inc.**

19-C Schoolhouse Road,  
Somerset, NJ 08873  
Phone: (732) 564.0700  
Fax: (732) 564.1515

U.S. Regional Office:  
4710B Interstate Drive,  
Cincinnati, Ohio 45264

Email: [USAsales@rosesci.com](mailto:USAsales@rosesci.com)

[www.rosescientific.com](http://www.rosescientific.com)