

# PVM Mixer



**PVM Shown can hold 2 Low or High Well Plates**

## **PULSING VORTEX MIXER**

- Variable speed: 100-2000 rpm
- Pulsing feature interrupts action for more thorough mixing
- Pulse rate: 60-80 ppm
- On/Off switch (power and pulse)
- Toggle clamp for securing deep well plates
- Contact us for custom clamping devices

### **Specifications:**

Electrical: 6' power cord (4' on 240V unit)  
Dimensions: 10" w x 7" d x 9" h, drive plate 7" x 7.75"  
Finish: Foxtip gray/charcoal  
Weight: 22 lbs.

Catalog Number	Price
099A PVM12 (120 volts)	\$1250.00
099A PVM24 (240 volts)	\$1270.00

### **Pulsing Feature:**

Glas-Col's pulsing feature is what makes these vortexers stand apart from other mixers. The pulsing action allows the vortex to collapse, before throwing back a full vortex. This time delayed sequence helps to join the upper and lower strata for more thorough and reproducible mixing.

REPEATABLE.

RELIABLE.

RESULTS.

### **Glas-Col, LLC**

711 Hulman Street., P.O. Box 2128, Terre Haute, IN 47802-0128 USA

Phone. (812) 235-6167, Fax. (812) 234-6975

E-mail: [pinnacle@glascol.com](mailto:pinnacle@glascol.com)

[www.glascol.com](http://www.glascol.com)



## **Rose Scientific Inc.**

19-C Schoolhouse Road,  
Somerset, NJ 08873  
Phone: (732) 564.0700  
Fax: (732) 564.1515

U.S. Regional Office:  
4710B Interstate Drive,  
Cincinnati, Ohio 45264

Email: [USAsales@rosesci.com](mailto:USAsales@rosesci.com)

[www.rosescientific.com](http://www.rosescientific.com)