

Catalog Number	Description	Product Weight (lbs)
099A WS18012	180° shaker with two finger clamps, 120V	11
099A WS18024	180° shaker with two finger clamps, 240V	11

WS 180° SHAKER

This shaker efficiently agitates small to medium solid phase peptide synthesis flasks. It is also useful for shaking culture flasks and for agitating hard-to-dissolve solids such as polymers.

- Shakes two flasks simultaneously.
- Simulates shaking by hand by rotating through 180° of arc.
- May be used with screw-cap top or side-arm flasks.
- Two 3-finger clamps included to secure glassware up to 12mm in shaft diameter.
- Shaft bearing design permits operation in cold room or incubator.
- DC filtering promotes longer motor life.
- Base configuration prevents "walking" on bench top

Specifications

Electrical: Fuse protector; on/off switch; variable speed control; 6-ft., 3-wire grounded cord.

Electrical Rating: 120 or 240V, 50/60HZ.

Speed Range: 2-30 rpm.

Size: 14"w x 6.5"d x 7"h.

Finish: Foxtip gray and charcoal.



19-C Schoolhouse Road,
Somerset, NJ 08873
Phone: (732) 564.0700
Fax: (732) 564.1515

U.S. Regional Office:
4710B Interstate Drive,
Cincinnati, Ohio 45264

Email: USAsales@rosesci.com

www.rosescientific.com